

# mpxtractor: A flexible R package to process data from different plate reader machines



**Student: Martin Banchero**



**Supervisor: Douwe Molenaar, PhD**

# Introduction

## High-throughput experiments

### Broad popularity:

- Drug discovery
  - Food research
  - Biofuels
- 
- Test hundreds of samples simultaneously.
  - Allow to collect, analyze data about complex interactions within biological systems.


96 well microplate






384 well microplate



dimensions of, 1536,  
3456 and more



- i.e Growth rates



- i.e Enzyme rates

# Introduction

## Microplate reader machines in the lab

**FluoStar Omega**



- Fluorescence intensity,
- Time-resolved fluorescence
- Absorbance
- Luminescence.

**Multiskan GO**



- Spectral scanning
- Single point in time (endpoint)
- Over specified period of time (Kinetic)

**SpectraMax**



- Optical density (OD)
- endpoint
- kinetic
- spectral scan



# Introduction

## Microplate reader machines: Raw data


### Output file FLUOStar Omega

- Measures at 6 time points

Schematic representation 96 microplate

```

1 User: USER      Path: C:\Program Files (x86)\BMG\Omega\User\Data\
2 Test ID
221
3 Test Name: Test_Fluorescence   Date: 29/08/2019      Time:
15:04:05
4 ID1: Reading
1
5 ID2:
Fluorescence
6 Fluorescence
(FI)
7
8
9 "Well
10 Row"      "Well
11 Col"      Content      "Raw Data (485, 520)
12 1 - 0 h   "      "Raw Data (485, 520)
13 2 - 0 h 2 min "      "Raw Data (485, 520)
14 3 - 0 h 7 min "      "Raw Data (485, 520)
15 4 - 0 h 12 min "      "Raw Data (485, 520)
16 5 - 0 h 17 min "      "Raw Data (485, 520)
17 6 - 0 h 22 min "      "Raw Data (485, 520)
18          0          0.043333333      0.125555556      0.21      0.292222222      0.376666667
19 A      1      Sample X1      186705      186950      188978      189018      189620      191198
20 A      2      Sample X2      186784      186224      187730      187383      188480      188519
  
```



Fluorescence



# Introduction

## Microplate reader machines: Raw data

### Output file Multiskan GO: List format

1	Plate	Well	Group	Type	Sample	Wavelength	Reading	Abs	Meas. Time [s]
2	Plate 1	C02	Assay	Calibrator	1 1/1	340	1	0.231	0.330
3	Plate 1	C02	Assay	Calibrator	1 1/1	340	2	0.233	30.320
4	Plate 1	C02	Assay	Calibrator	1 1/1	340	3	0.235	60.320
5	Plate 1	C02	Assay	Calibrator	1 1/1	340	4	0.237	90.320
6	Plate 1	D02	Assay	Calibrator	2 1/1	340	1	0.236	0.690
7	Plate 1	D02	Assay	Calibrator	2 1/1	340	2	0.239	30.690
8	Plate 1	D02	Assay	Calibrator	2 1/1	340	3	0.240	60.690
9	Plate 1	D02	Assay	Calibrator	2 1/1	340	4	0.242	90.690
10	Plate 1	E02	Assay	Calibrator	3 1/1	340	1	0.241	1.020
11	Plate 1	E02	Assay	Calibrator	3 1/1	340	2	0.244	31.020
12	Plate 1	E02	Assay	Calibrator	3 1/1	340	3	0.246	61.020
13	Plate 1	E02	Assay	Calibrator	3 1/1	340	4	0.250	91.020
14	Plate 1	F02	Assay	Calibrator	4 1/1	340	1	0.231	1.390
15	Plate 1	F02	Assay	Calibrator	4 1/1	340	2	0.235	31.380
16	Plate 1	F02	Assay	Calibrator	4 1/1	340	3	0.238	61.380
17	Plate 1	F02	Assay	Calibrator	4 1/1	340	4	0.242	91.380
18	Plate 1	C03	Assay	Control	5 1/1	340	1	0.240	0.350
19	Plate 1	C03	Assay	Control	5 1/1	340	2	0.246	30.340
20	Plate 1	C03	Assay	Control	5 1/1	340	3	0.251	60.340
21	Plate 1	C03	Assay	Control	5 1/1	340	4	0.257	90.340
22	Plate 1	D03	Assay	Control	6 1/1	340	1	0.240	0.680

### Output file Multiskan GO: Table format

1	Results of Photometric1																			
2	Plate: Plate 1 - Wavelength: 600																			
3	Reading: 1																			
4	0.081	0.081	0.081	0.080	0.080	0.082	0.081	0.083	0.081	0.083	0.081	0.082	0.084	0.082	0.082	0.081	0.082	0.081	0.082	0.082
5	0.081	0.080	0.081	0.081	0.081	0.082	0.082	0.083	0.083	0.083	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
6	0.082	0.086	0.087	0.085	0.085	0.085	0.083	0.083	0.083	0.083	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
7	0.084	0.082	0.083	0.082	0.083	0.083	0.082	0.083	0.083	0.083	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
8	0.081	0.085	0.087	0.085	0.085	0.085	0.083	0.083	0.083	0.083	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
9	0.082	0.082	0.083	0.083	0.083	0.083	0.083	0.084	0.082	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
10	0.083	0.082	0.082	0.082	0.083	0.084	0.083	0.083	0.082	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
11	0.083	0.086	0.085	0.085	0.085	0.085	0.084	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
12	0.083	0.082	0.082	0.082	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
13	0.083	0.082	0.082	0.082	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
14	0.083	0.086	0.085	0.085	0.085	0.085	0.084	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
15	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
16	0.083	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
17	0.083	0.082	0.082	0.082	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
18	0.084	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
19	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
20	0.084	0.085	0.084	0.085	0.085	0.085	0.084	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
21	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
22	0.086	0.084	0.084	0.086	0.086	0.086	0.087	0.086	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085
23	0.083	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
24	0.083	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085	0.085
25	0.086	0.084	0.083	0.085	0.087	0.086	0.087	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
26	0.082	0.086	0.086	0.087	0.087	0.086	0.085	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
27	0.085	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
28	0.082	0.085	0.086	0.087	0.087	0.086	0.085	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
29	0.085	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084	0.084
30	0.082	0.085	0.086	0.087	0.087	0.086	0.085	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086	0.086
31	0.085	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
32	0.084	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083	0.083
33	0.079	0.080	0.079	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081
34	0.080	0.080	0.080	0.082	0.082	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081	0.081

● 4 readings every 30 seconds

● No time present in table format

1	Results of Photometric1									
2	Plate: Plate 1 - Wavelength: 600									
3	Reading: 1									
4	0.088	0.090	0.090	0.090	0.088	0.090	0.091	0.090	0.090	0.089
5	0.091	0.093	0.090	0.090	0.089	0.089	0.089	0.088	0.088	0.089
6	0.088	0.089	0.089							



# Introduction

## Microplate reader machines: Raw data

### Output file SpectraMax(.txt)

Time(hh:mm:ss)	Temperature(°C)	A1	A2	A3	A4	A5	A6	A7	A8		
A17	A18	A19	A20	A21	A22	A23	A24	B1	B2	B3	B4
B13	B14	B15	B16	B17	B18	B19	B20	B21	B22	B23	B24
C9	C10	C11	C12	C13	C14	C15	C16	C17	C18	C19	C20
D5	D6	D7	D8	D9	D10	D11	D12	D13	D14	D15	D16
E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12
E21	E22	E23	E24	F1	F2	F3	F4	F5	F6	F7	F8
F17	F18	F19	F20	F21	F22	F23	F24	G1	G2	G3	G4
G13	G14	G15	G16	G17	G18	G19	G20	G21	G22	G23	G24
H9	H10	H11	H12	H13	H14	H15	H16	H17	H18	H19	H20
I5	I6	I7	I8	I9	I10	I11	I12	I13	I14	I15	I16
J1	J2	J3	J4	J5	J6	J7	J8	J9	J10	J11	J12
J21	J22	J23	J24	K1	K2	K3	K4	K5	K6	K7	K8
K17	K18	K19	K20	K21	K22	K23	K24	L1	L2	L3	L4
L13	L14	L15	L16	L17	L18	L19	L20	L21	L22	L23	L24
M9	M10	M11	M12	M13	M14	M15	M16	M17	M18	M19	M20
N5	N6	N7	N8	N9	N10	N11	N12	N13	N14	N15	N16
O1	O2	O3	O4	O5	O6	O7	O8	O9	O10	O11	O12
O21	O22	O23	O24	P1	P2	P3	P4	P5	P6	P7	P8
P17	P18	P19	P20	P21	P22	P23	P24				
0:00	30,00	0,0821	0,1387	0,1438	0,1558	0,1416	0,1425	0,15	0,144	0,1404	0,1396
0,1526	0,1343	0,1551	0,1399	0,1458	0,0833	0,1552	0,2537	0,2535	0,2781	0,2708	0,2755
0,2769	0,2466	0,2459	0,2505	0,2529	0,2623	0,2626	0,2708	0,2559	0,134	0,1536	0,2529
0,3049	0,2984	0,2645	0,2621	0,26	0,2732	0,2662	0,2632	0,2606	0,2302	0,2265	0,2956
0,2775	0,2796	0,2774	0,2822	0,2748	0,2748	0,2737	0,2878	0,2792	0,2786	0,273	0,282
0,1655	0,17	0,1709	0,1728	0,1718	0,1803	0,1702	0,171	0,1694	0,1779	0,1707	0,1744
0,1715	0,1757	0,1418	0,2444	0,2708	0,2746	0,2765	0,2851	0,2832	0,2881	0,2845	0,284
0,2801	0,2912	0,294	0,296	0,2817	0,1929	0,1425	0,2675	0,2616	0,2755	0,281	0,2843

- Different researchers using different formats, .txt, .csv, xlsx

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE	AF	AG	AH	AI	
1	Time(hh:mm:ss)	Temperature(°C)	A1	A2	A3	A4	A5	A6	A7	A8	A9				B1	B2	B3	B4	B5	B6	B7	B8	B9				B12	C1	C2	C3	C4	C5	C6	C7	C8	C9
2	0:00	37							0.3911	0.3963	0.3995										0.7936	0.6018	0.8342										1.2167	0.8148	1.2405	
3	0:05	37							0.3902	0.3969	0.4003										0.7962	0.6028	0.836										1.2196	0.8231	1.2427	
4	0:10	37							0.3908	0.3969	0.4011										0.7971	0.6031	0.8368										1.221	0.8243	1.2438	
5	0:15	37							0.3907	0.3969	0.4013										0.7979	0.6034	0.8374										1.2218	0.8259	1.2442	
6	0:20	37							0.3909	0.3969	0.4007										0.7985	0.6039	0.8377										1.2234	0.8284	1.2459	
7	0:25	37							0.3909	0.3973	0.4004										0.7993	0.6042	0.8381										1.2241	0.8289	1.2464	
8	0:30	37							0.391	0.3974	0.4008										0.8	0.6046	0.8388										1.2252	0.8299	1.247	
9	0:35	37							0.3911	0.3975	0.4008										0.8004	0.6048	0.84										1.2255	0.8307	1.2471	

# Introduction

## Experimental design: Layout files

## Layout files contain:

- A description of the location and content of each well.
- Indicating if the well is control or a sample
- Control variables, like concentration or pH.
- Information in layout files not present in raw output files from plate readers.

[illegible]

# Introduction

## Recap:

- **Different possible output files formats.**
  - .txt, csv, xlsx depends on the user
  - Different data representation, i.e multiskan GO list or table format depending on the manufacturer.
- **Information in layout files not included in raw data files.**



# Introduction

## Previous packages to process microplate raw data

- **cellHTS2(2006)**

- Framework for RNAi screening
- Raw files has to be formatted before analysis
- *Ref: M. Boutros, L. P. Brás, and W. Huber, “Analysis of cell-based RNAi screens,” Genome Biol.*

- **Douwe first package(2010)**

- Parser only for SpectraMax, also need for update

- **Plater package(2016)**

- Designed to read plate-shaped file format of experimental design into tidy data frame.
- Also, allows the conversion of a tidy data frame into plate-shaped data.
- *Ref: S. M. Hughes, “plater: Read, Tidy, and Display Data from Microtiter Plates.”*

# Introduction

## Tidy data

- Tidy data is a standard way to organize data
- Well suited for vectorised programming languages like R
- Easy to plug into tidyverse collection of R packages for data analysis

### Variables

```
#> # A tibble: 6 x 4
#>   Wells Time Temperature Measurement
#>   <chr> <chr> <dbl> <dbl>
#> 1 A01 00:00:00 30,00 0.0821
#> 2 A01 00:02:00 30,00 0.0811
#> 3 A01 00:04:00 30,00 0.081
#> 4 A02 00:00:00 30,00 0.089
#> 5 A02 00:02:00 30,00 0.0818
#> 6 A02 00:04:00 30,00 0.0822
```

- Each variable represented as column

### Observations

```
#> # A tibble: 6 x 4
#>   Wells Time Temperature Measurement
#>   <chr> <chr> <dbl> <dbl>
#> 1 A01 00:00:00 30,00 0.0821
#> 2 A01 00:02:00 30,00 0.0811
#> 3 A01 00:04:00 30,00 0.081
#> 4 A02 00:00:00 30,00 0.089
#> 5 A02 00:02:00 30,00 0.0818
#> 6 A02 00:04:00 30,00 0.0822
```

- Each observation represented by a row

### Values

```
#> # A tibble: 6 x 4
#>   Wells Time Temperature Measurement
#>   <chr> <chr> <dbl> <dbl>
#> 1 A01 00:00:00 30,00 0.0821
#> 2 A01 00:02:00 30,00 0.0811
#> 3 A01 00:04:00 30,00 0.081
#> 4 A02 00:00:00 30,00 0.089
#> 5 A02 00:02:00 30,00 0.0818
#> 6 A02 00:04:00 30,00 0.0822
```

- Each value belong to one cell

# Introduction

## Project goal

**The main goal of this research project is to develop an R package to:**

- Reduce the time of data pre-processing in high-throughput experiments and improve standardization.
  - Process different output files from micro-plates readers. Store this information as tidy data.
  - Process layout files. Store this into tidy data
  - Combine output files with layout files into one tidy tibble object.

# Introduction

## mpxtractor

- The main function of mpxtractor is to provide a simple method that can process raw data from microplate readers:
  - SpectraMax
  - FluorStar
  - MultiscanGO
- This package generates tidy tibble objects, this data can be easily combined with layout files.
- Growth rates can be calculated and plot over a microplate-shape format.



Package available here:

<https://github.com/MartinBanchero/mpxtractor>

**open source license**  
**GPL-3**



# Data

## **SpectraMax**

2 files:

Time: 26 hs and  
16 hs, interval 2  
minutes. Temperature:  
30°C

## **Multiskan GO:**

2 files:

1000 readings  
Time interval 2  
minutes.  
Demo files for list  
format.

## **FluorStar:**

2 files:

Time 20 hs 30 minutes  
Time interval 2 minutes




Sieze Douwenga, PhD student at  
Systems Biology lab-VU

# Methods

## Mpxtractor: flowchart

- Wrangling functions
- Plotting functions




# Methods

## Wrangling functions: SpectraMax files

### ● Caption raw file SpectraMax

Time(hh:mm:ss)	Temperature(°C)	A1	A2	A3	A4	A5	A6	A7	A8
A17	A18	A19	A20	A21	A22	A23	A24	B1	B2
B13	B14	B15	B16	B17	B18	B19	B20	B21	B22
C9	C10	C11	C12	C13	C14	C15	C16	C17	C18
D5	D6	D7	D8	D9	D10	D11	D12	D13	D14
E1	E2	E3	E4	E5	E6	E7	E8	E9	E10
E21	E22	E23	E24	F1	F2	F3	F4	F5	F6
F17	F18	F19	F20	F21	F22	F23	F24	G1	G2
G13	G14	G15	G16	G17	G18	G19	G20	G21	G22
H9	H10	H11	H12	H13	H14	H15	H16	H17	H18
I5	I6	I7	I8	I9	I10	I11	I12	I13	I14
J1	J2	J3	J4	J5	J6	J7	J8	J9	J10
J21	J22	J23	J24	K1	K2	K3	K4	K5	K6
K17	K18	K19	K20	K21	K22	K23	K24	L1	L2
L13	L14	L15	L16	L17	L18	L19	L20	L21	L22
M9	M10	M11	M12	M13	M14	M15	M16	M17	M18
N5	N6	N7	N8	N9	N10	N11	N12	N13	N14
O1	O2	O3	O4	O5	O6	O7	O8	O9	O10
O21	O22	O23	O24	P1	P2	P3	P4	P5	P6
P17	P18	P19	P20	P21	P22	P23	P24		
00:00	30,00	0,0821	0,1387	0,1438	0,1558	0,1416	0,1425	0,15	0,144
0,1536	0,1343	0,1399	0,1399	0,1458	0,0833	0,1552	0,2537	0,2535	0,2781
0,2769	0,2466	0,2459	0,2585	0,2529	0,2623	0,2626	0,2708	0,2559	0,134
0,3049	0,2984	0,2645	0,2621	0,26	0,2732	0,2662	0,2632	0,2606	0,2302
0,2775	0,2796	0,2774	0,2822	0,2748	0,2748	0,2737	0,2878	0,2792	0,2786
0,1655	0,17	0,1709	0,1728	0,1718	0,1803	0,1702	0,171	0,1694	0,1779
0,1715	0,1757	0,1418	0,2444	0,2708	0,2746	0,2765	0,2851	0,2832	0,2881
0,2801	0,2912	0,294	0,296	0,2817	0,1020	0,1425	0,2616	0,2755	0,281



### ● Example of usage and output of read\_spectramax\_data()

```
R> data_spectramax <- mpxtractor::read_spectramax_data(file = file_path)
R> head(data_spectramax)
```

```
## # A tibble: 6 x 4
##   Wells Time Temperature Measurement
##   <chr> <chr> <chr> <dbl>
## 1 A01 00:00:00 30,00 0.0821
## 2 A01 00:02:00 30,00 0.0811
## 3 A01 00:04:00 30,00 0.081
## 4 A02 00:00:00 30,00 0.139
## 5 A02 00:02:00 30,00 0.0818
## 6 A02 00:04:00 30,00 0.0822
```

# Methods

## Wrangling functions: Multiscan GO files

### ● Caption raw file: table format

1 Results of Photometric1

2 Plate: Plate 1 - Wavelength: 600

3

4

5 Reading: 1

0.081	0.081	0.080	0.080	0.082	0.081	0.083	0.081	0.083	0.081	0.082	0.084	0.082	0.082	0.081	0.082	0.081
0.082	0.086	0.087	0.085	0.085	0.083	0.083	0.083	0.084	0.084	0.082	0.086	0.082	0.084	0.084	0.085	0.084
0.081	0.085	0.087	0.085	0.085	0.083	0.083	0.083	0.083	0.083	0.082	0.086	0.084	0.082	0.083	0.084	0.082
0.082	0.086	0.085	0.085	0.082	0.082	0.083	0.083	0.083	0.083	0.082	0.085	0.082	0.082	0.083	0.082	0.083
0.083	0.086	0.085	0.084	0.083	0.084	0.083	0.084	0.084	0.083	0.082	0.086	0.083	0.083	0.083	0.083	0.083
0.083	0.085	0.084	0.085	0.084	0.084	0.083	0.083	0.083	0.083	0.085	0.083	0.082	0.083	0.083	0.082	0.082
0.084	0.086	0.086	0.085	0.085	0.085	0.084	0.085	0.085	0.084	0.083	0.086	0.087	0.083	0.083	0.084	0.083
0.084	0.085	0.084	0.085	0.085	0.084	0.083	0.085	0.084	0.083	0.085	0.084	0.083	0.083	0.083	0.082	0.083
0.084	0.087	0.085	0.087	0.086	0.086	0.085	0.085	0.086	0.086	0.086	0.087	0.086	0.083	0.084	0.085	0.086
0.083	0.086	0.086	0.086	0.087	0.085	0.085	0.085	0.084	0.084	0.086	0.086	0.083	0.085	0.083	0.084	0.085
0.084	0.086	0.085	0.085	0.087	0.084	0.086	0.085	0.085	0.086	0.086	0.085	0.087	0.086	0.084	0.085	0.086
0.083	0.084	0.087	0.086	0.087	0.086	0.083	0.086	0.086	0.086	0.086	0.086	0.086	0.085	0.085	0.084	0.086
0.082	0.086	0.086	0.087	0.086	0.085	0.086	0.086	0.087	0.087	0.086	0.085	0.086	0.083	0.083	0.086	0.086
0.082	0.085	0.086	0.087	0.087	0.086	0.083	0.086	0.084	0.085	0.085	0.085	0.084	0.083	0.083	0.083	0.086
0.082	0.084	0.084	0.085	0.085	0.085	0.084	0.083	0.084	0.083	0.085	0.087	0.084	0.082	0.084	0.083	0.084
0.079	0.080	0.079	0.081	0.081	0.081	0.079	0.080	0.080	0.080	0.081	0.081	0.081	0.081	0.081	0.080	0.080

22


23 Reading: 2

0.081	0.082	0.081	0.080	0.082	0.081	0.082	0.081	0.083	0.081	0.082	0.084	0.082	0.082	0.081	0.082	0.081
0.082	0.086	0.087	0.085	0.085	0.084	0.083	0.084	0.084	0.084	0.083	0.086	0.083	0.084	0.085	0.085	0.084

```
R> data_multiscan <- mpctractor::read_multiscango_data(file = file_path,  
time_interval = "2.5 min")  
R> head(data_multiscan)
```

```
#> # A tibble: 6 x 5  
#> Wells Time Reading wavelength Measurement  
#> <chr> <chr> <dbl> <dbl> <dbl>  
#> 1 A01 00:00:00 1 600 0.081  
#> 2 A01 00:02:30 2 600 0.081  
#> 3 A01 00:05:00 3 600 0.081  
#> 4 A01 00:07:30 4 600 0.081  
#> 5 A01 00:10:00 5 600 0.081  
#> 6 A01 00:12:30 6 600 0.081
```

### ● Example of usage and output of read\_multiscan\_data()






# Methods

## Wrangling functions: Multiscan GO files

### ● Caption raw file: List format

	Plate	Well	Group	Type	Sample	Wavelength	Reading	Abs	Meas. Time [s]
1	Plate 1	C02	Assay	Calibrator	1 1/1	340	1	0.231	0.330
2	Plate 1	C02	Assay	Calibrator	1 1/1	340	2	0.231	30.320
3	Plate 1	C02	Assay	Calibrator	1 1/1	340	3	0.231	60.320
4	Plate 1	C02	Assay	Calibrator	1 1/1	340	4	0.231	90.320
5	Plate 1	D02	Assay	Calibrator	2 1/1	340	1	0.236	0.690
6	Plate 1	D02	Assay	Calibrator	2 1/1	340	2	0.236	30.690
7	Plate 1	D02	Assay	Calibrator	2 1/1	340	3	0.240	60.690
8	Plate 1	D02	Assay	Calibrator	2 1/1	340	4	0.240	90.690
9	Plate 1	E02	Assay	Calibrator	3 1/1	340	1	0.241	1.020
10	Plate 1	E02	Assay	Calibrator	3 1/1	340	2	0.244	31.020
11	Plate 1	E02	Assay	Calibrator	3 1/1	340	3	0.246	61.020
12	Plate 1	E02	Assay	Calibrator	3 1/1	340	4	0.250	91.020
13	Plate 1	F02	Assay	Calibrator	4 1/1	340	1	0.231	1.390
14	Plate 1	F02	Assay	Calibrator	4 1/1	340	2	0.235	31.380
15	Plate 1	F02	Assay	Calibrator	4 1/1	340	3	0.238	61.380
16	Plate 1	F02	Assay	Calibrator	4 1/1	340	4	0.242	91.380
17	Plate 1	C03	Assay	Control	5 1/1	340	1	0.240	0.350
18	Plate 1	C03	Assay	Control	5 1/1	340	2	0.246	30.340
19	Plate 1	C03	Assay	Control	5 1/1	340	3	0.251	60.340
20	Plate 1	C03	Assay	Control	5 1/1	340	4	0.257	90.340
21	Plate 1	D03	Assay	Control	6 1/1	340	1	0.240	0.600
22	Plate 1	D03	Assay	Control	6 1/1	340	2	0.246	30.670
23	Plate 1	D03	Assay	Control	6 1/1	340	3	0.252	60.670
24	Plate 1	D03	Assay	Control	6 1/1	340	4	0.259	90.670
25	Plate 1	E03	Assay	Control	7 1/1	340	1	0.274	1.040
26	Plate 1	E03	Assay	Control	7 1/1	340	2	0.281	31.030
27	Plate 1	E03	Assay	Control	7 1/1	340	3	0.289	61.030
28	Plate 1	E03	Assay	Control	7 1/1	340	4	0.296	91.040
29	Plate 1	B02	Assay	Blank	Blank_Assay 1/1	340	1	0.231	0.000
30	Plate 1	B02	Assay	Blank	Blank_Assay 1/1	340	2	0.230	30.000
31	Plate 1	B02	Assay	Blank	Blank_Assay 1/1	340	3	0.231	60.000
32	Plate 1	B02	Assay	Blank	Blank_Assay 1/1	340	4	0.231	90.000
33	Plate 1	B02	Assay	Blank	Blank_Assay 1/1	340	4	0.231	90.000



```
R> data_multiscan <- mpctractor::read_multiscango_data(file_path, input_type = "list")
R> head(data_multiscan)
```

Wells	Time	Measurement	Plate	Group	Type	Sample	Wavelength	Reading	Meas.time.sec
<chr>	<chr>	<dbl>	<chr>	<chr>	<chr>	<chr>	<dbl>	<dbl>	<dbl>
C02	00:00:00	0.231	Plate 1	Assay	Calibrator 1	1/1	340	1	0.33
C02	00:00:30	0.233	Plate 1	Assay	Calibrator 1	1/1	340	2	30.32
C02	00:01:00	0.235	Plate 1	Assay	Calibrator 1	1/1	340	3	60.32
C02	00:01:30	0.237	Plate 1	Assay	Calibrator 1	1/1	340	4	90.32
D02	00:00:00	0.236	Plate 1	Assay	Calibrator 2	1/1	340	1	0.69
D02	00:00:30	0.239	Plate 1	Assay	Calibrator 2	1/1	340	2	30.69

### ● Example of usage and output of read\_multiscan\_data()

# Methods

## Wrangling functions: FluorStar files

- **Caption raw file**

```


User: USER                               Path: C:\Program Files (x86)\BNC\Omega\User\Data\
221
2 Test Name: Test_Fluorescence             Date: 29/08/2019             Time:
15:04:05
3 ID1: Reading
1
4 ID2:
Fluorescence
5 Fluorescence
(#1)
6
7
8 Well
9 Row
Content
"Raw Data (485, 520)"
10 Col
11 1 0 h 1 min "Raw Data (485, 520)"
12 2 0 h 2 min "Raw Data (485, 520)"
13 3 0 h 7 min "Raw Data (485, 520)"
14 4 0 h 12 min "Raw Data (485, 520)"
15 5 0 h 17 min "Raw Data (485, 520)"
16 6 0 h 22 min "Raw Data (485, 520)"
17
18 A 1 Sample X1 0.04333333 0.12555556 0.21 0.29222222 0.37666667
186705 186950 188978 189018 189620 191198

```

## ● Example of usage and output of `read_fluorstar_data()`

```
R> data_fluorstar <- mp_extractor::read_fluorstar_data(file = file_path)
R> head(data_fluorstar)
```

```
#> # A tibble: 6 x 5
#>   Wells Sample      Time Measurement wavelength
#>   <chr> <chr>    <chr>    <int>      <chr>
#> 1 A01 " X1" 00:00:00 186705 485, 520
#> 2 A01 " X1" 00:02:00 186950 485, 520
#> 3 A01 " X1" 00:07:00 188978 485, 520
#> 4 A01 " X1" 00:12:00 189018 485, 520
#> 5 A01 " X1" 00:17:00 189620 485, 520
#> 6 A01 " X1" 00:22:00 191198 485, 520
```



# Methods

## Read multiple data files


### Usage and output of read\_multiple\_data\_files

```
R> data_multiple_spectramax <- mpextractor::read_multiple_data_files(  
  reader_type = "spectramax",  
  filename = file_path,  
  plate_names = c("plate_1", "plate_2")  
)  
R> head(data_multiple_spectramax)
```

```
#> # A tibble: 6 x 5  
#>   Wells Time Temperature Measurement plate_filename  
#>   <chr> <chr>      <dbl>      <dbl>      <chr>  
#> 1 A01 00:00:00 30,00      0.0821    plate_1  
#> 2 A01 00:02:00 30,00      0.0811    plate_1  
#> 3 A01 00:04:00 30,00      0.081     plate_1  
#> 4 A02 00:00:00 30,00      0.139     plate_1  
#> 5 A02 00:02:00 30,00      0.0818    plate_1  
#> 6 A02 00:04:00 30,00      0.0822    plate_1
```

```
R> tail(data_multiple_spectramax)
```

```
#> # A tibble: 6 x 5  
#>   Wells Time Temperature Measurement plate_filename  
#>   <chr> <chr>      <dbl>      <dbl>      <chr>  
#> 1 P23 00:00:00 30,00      0.0801    plate_2  
#> 2 P23 00:02:00 30,00      0.0799    plate_2  
#> 3 P23 00:04:00 30,00      0.08      plate_2  
#> 4 P24 00:00:00 30,00      0.0812    plate_2  
#> 5 P24 00:02:00 30,00      0.0811    plate_2  
#> 6 P24 00:04:00 30,00      0.0815    plate_2
```



This function takes as argument the type of reader machine (blue) and the name of the file (red)

# Methods

## Read layout files

### ● Caption layout file


	A	B	C	D	E	F	G	H	I	J	K	L	M
1	basic	1	2	3	4	5	6	7	8	9	10	11	12
2	A	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
3	B	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
4	C	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
5	D	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
6	E	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
7	F	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
8	G	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
9	H	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	control
11	medium	1	2	3	4	5	6	7	8	9	10	11	12
12	A	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
13	B	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
14	C	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
15	D	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
16	E	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
17	F	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
18	G	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
19	H	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup	fit_milk_sup
21	concentration	1	2	3	4	5	6	7	8	9	10	11	12
22	A	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
23	B	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
24	C	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
25	D	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
26	E	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
27	F	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
28	G	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0
29	H	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12E-05	1.024E-05	0

```
R> df_layout <- mpxtractor::read_layout_file(file)
```

```
R> head(df_layout)
```

	Wells	basic	medium	concentration
	<chr>	<chr>	<chr>	<dbl>
1	A01	control	milk	100
2	A02	sample	milk	20
3	A03	sample	milk	4
4	A04	sample	milk	0.8
5	A05	sample	milk_sup	0.16
6	A06	sample	milk_sup	0.032

### ● Example of usage and output of read\_layout\_file()



**read\_layout\_file()** similar to  
**read\_plate()** from plater package

S. M. Hughes, “*plater: Read, Tidy, and Display Data from Microtiter Plates.*”

# Methods

## Combine data files with layout files


**df\_data:** Raw data into tibble

**Reader\_type:** Specified reader machine

```
R> df_data_fl <- mpextractor::read_fluorstar_data(file_fl)
R> df_data_layout <- mpextractor::combine_data_with_layout(
  df_data = df_data_fl,
  reader_type = "fluorstar",
  layout_files = file
)
R> head(df_data_layout)
```

A tibble: 6 x 9								
Wells	Sample	Time	Measurement	wavelength	basic	medium	concentration	layout_file
<chr>	<chr>	<chr>	<dbl>	<chr>	<chr>	<chr>	<dbl>	<chr>
1 A01	* X1*	00:00:00	3.5	600	control	milk	100	test_fluorstar_layout
2 A01	* X1*	00:02:00	3.44	600	control	milk	100	test_fluorstar_layout
3 A01	* X1*	00:07:00	3.5	600	control	milk	100	test_fluorstar_layout
4 A01	* X1*	00:12:00	3.5	600	control	milk	100	test_fluorstar_layout
5 A01	* X1*	00:17:00	3.21	600	control	milk	100	test_fluorstar_layout
6 A01	* X1*	00:22:00	3.35	600	control	milk	100	test_fluorstar_layout


**Raw data combined with layout file**





# Methods

## Plot layout file



medium

● filt\_milk\_sup


● milk

● milk\_sup

basic

● control

▲ sample





## Layout file

	A	B	C	D	E	F	G	H	I	J	K	L	M
1	basic	1	2	3	4	5	6	7	8	9	10	11	12
2	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
3	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
4	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
5	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
6	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
7	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
8	control	sample	sample	sample	sample	sample	sample	control	sample	sample	sample	sample	control
9	medium	1	2	3	4	5	6	7	8	9	10	11	12
10	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
11	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
12	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
13	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
14	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
15	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
16	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
17	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
18	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
19	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
20	milk	milk	milk	milk	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup	milk_sup
21	concentration	1	2	3	4	5	6	7	8	9	10	11	12
22	A	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
23	B	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
24	C	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
25	D	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
26	E	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
27	F	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
28	G	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0
29	H	100	20	4	0.8	0.16	0.032	0.0064	0.00128	0.000256	5.12e-05	1.024e-05	0

Representation up to 2 variables through shape, color and concentration(optional)

# Methods

## Plot growth rates



```
R> growth_rate_plot<- mpxtractor::plot_gr_microplate(
df_data = df_corrected,
var_gr = "Measurement",
cond_to_col = "condition",
plate_file = NULL,
ws = "2hs",
output_filename = NULL)
```

# Methods


## Plot growth rates

- **compute\_growth\_rates()**

This function calculate growth rate using **Savitzky and Golay filter** from package Signal.

- **Sav-Gol filter :**

- Type of low pass filter used to smooth the data, improve graphical output due to large number of data points and random noise.
- It can be used to compute smoothed first order derivatives to obtain the growth rates.
- Use a moving window average to fit a polynomial function through the data points





# Methods

## Plot growth rates

### Time series

```
#> # A tibble: 6 x 9
#>   Wells Sample Time Measurement wavelength Basic condition concentration layout_file
#>   <chr> <chr> <chr> <dbl>         <chr>      <chr> <chr>      <dbl>      <chr>
#> 1 A01 " X1" 00:00:00 3.5           600      poor    milk        100    test_fluorstar_layout
#> 2 A01 " X1" 00:02:00 3.44          600      poor    milk        100    test_fluorstar_layout
#> 3 A01 " X1" 00:07:00 3.5           600      poor    milk        100    test_fluorstar_layout
#> 4 A01 " X1" 00:12:00 3.5           600      poor    milk        100    test_fluorstar_layout
#> 5 A01 " X1" 00:17:00 3.21          600      poor    milk        100    test_fluorstar_layout
#> 6 A01 " X1" 00:22:00 3.35          600      poor    milk        100    test_fluorstar_layout
```

Time series has to be fixed by user before calculate the growth rates, otherwise mpctractor through an error.

### Remark: missing values

- Missing values in column “Measurements” are imputed by taking the average between surrounding numbers.
- In this case mpctractor shows a warning letting the user know about the imputation.

# Results

## Validation: Growth rates

### Fluorstar data

Wells	mpxtractor data			Validation data		
	Time(hs )	log(Abs)	Growth rate	Time(h s)	log(Abs)	Growth rate
A01	0.03333	-0.84373	-0.07818094	0.00000	-0.843738	-0.078181
:	:	:	:	:	:	:
A01	1.28333	-1.175061	0.1089509	1.25000	-1.175061	0.108951
A01	1.36666	-0.717644	0.06203398	1.33333	-0.717645	-inf
:	:	:	:	:	:	:
A01	2.28333	-1.161871	-0.05384393	2.25000	-1.16187	-inf
A01	2.36666	-1.010418	0.057739	2.33333	-inf	-inf
A01	2.45000	-0.858965	0.00011882	2.41666	-0.858966	inf
:	:	:	:	:	:	:
A01	3.36666	-0.999400	0.06919943	3.33333	-0.999401	inf
A01	3.45000	-0.857078	0.08525832	3.41666	-0.857079	0.085258
:	:	:	:	:	:	:
A01	16.1166	-0.797175	-0.1556584	16.0833	-0.797175	0.073194
A01	--	--	--	20.4166	-0.920299	0.057739

Ws = 2hs

Value corresponding to background absorbance

Values imputed by taking the average of the surrounding values

# Results

## Validation: Growth rates

### Multiskan Go data

Wells	Time(hs)	mpxtractor data		Validation data	
		log(Abs)	Growth rate	log(Abs)	Growth rate
A01	0.000000	-6.907755	0.6355529	-Inf	NaN
:	:	:	:	:	:
A01	0.666667	-6.907755	0.6355529	-Inf	NaN
A01	0.708333	-6.907755	0.6355529	-6.907755	NaN
:	:	:	:	:	:
A01	0.958333	-6.907755	0.6355529	-6.907755	NaN
A01	1.000000	-6.907755	0.6355529	-6.907755	Inf
:	:	:	:	:	:
A01	1.666667	-6.214608	0.8386713	-6.214608	Inf
A01	1.708333	-5.809143	0.8284531	-5.809143	0.828453
:	:	:	:	:	:
A01	41.625000	-4.710531	-0.0173069	-4.710531	-0.0173069

Value corresponding to background absorbance

Values imputed by taking the average of the surrounding values

# Results

## Validation: Growth rates

### SpectraMax data

Wells	Time(hs)	mpxtractor data		Validation data	
		log(Abs)	Growth rate	log(Abs)	Growth rate
A01	0.0000	-7.264430	-0.238977	-7.264430	-0.238978
:	:	:	:	:	:
A01	1.1000	-6.907755	-0.384584	-6.907755	-0.384584
A01	1.1333	-6.907755	-0.412485	-6.907755	-Inf
:	:	:	:	:	:
A01	2.1000	-7.600902	0.033795	-7.600902	-Inf
A01	2.1333	-7.600902	0.031808	-Inf	-Inf
A01	2.1666	-7.600902	0.007807	-7.600902	Inf
:	:	:	:	:	:
A01	3.1333	-7.130899	0.357646	-7.130899	Inf
A01	3.1666	-7.600902	0.339011	-7.600902	0.339011
A01	26.1333	-6.645391	0.032598	-6.645391	0.032599

Value corresponding to background absorbance

Values imputed by taking the average of the surrounding values

# Conclusion

- **Reduction** of the time to process raw files into tidy data.
- Improvement in **standardization**, only two types of files.
- plate readers, acceptable number of different plate machines from three different
- taking the mean of the surrounding data points. Yes in the data and imputes by
- Help to improve **reproducibility**

# Future improvements

- Improve robustness of **mpxtractor** through users' reports of bugs in the package repository.
- More testing
- **Addition of features:**
  - To process data from other machines
  - Calculate enzyme rates
- Increase user friendliness using R-shiny

# Acknowledgments

## Systems Biology Lab-VU:



**Supervisor:** Douwe Molenaar PhD



Sieze Douwenga, PhD student



Vera Benavente, Lab technician

Dennis Botman, Phd student  
Laura Guilherme Luzia, Phd student

**Einde**



**Questions?**

# supplementary

## Set of parameters

Type of argument	Value	R parameter	Python parameter
Window length = window size / diff time	wl= (2/ 0.0333) hs	n=60	window_length=60
Filter order	1	p=1	polyorder=1
Derivative order	1	m=1	derivative=1

# supplementary

Input raw file for cellHTS with  
fluorescence values.

1	FT01-G01	A01	887763
2	FT01-G01	A02	958308
3	FT01-G01	A03	1012685
4	FT01-G01	A04	872603
5	FT01-G01	A05	1179875
6	FT01-G01	A06	1213279
7	FT01-G01	A07	1182103
8	FT01-G01	A08	1353041
9	FT01-G01	A09	1236055
10	FT01-G01	A10	1228895
11	FT01-G01	A11	1308768
12	FT01-G01	A12	1313472
13	FT01-G01	A13	1404374